

# Package: DiscoRhythm (via r-universe)

May 29, 2026

**Title** Interactive Workflow for Discovering Rhythmicity in Biological Data

**Version** 1.28.0

**Description** Set of functions for estimation of cyclical characteristics, such as period, phase, amplitude, and statistical significance in large temporal datasets. Supporting functions are available for quality control, dimensionality reduction, spectral analysis, and analysis of experimental replicates. Contains a R Shiny web interface to execute all workflow steps.

**Depends** R (>= 3.6.0)

**License** GPL-3

**Encoding** UTF-8

**LazyData** true

**RoxygenNote** 7.1.1

**VignetteBuilder** knitr

**Imports** matrixTests, matrixStats, MetaCycle (>= 1.2.0), data.table, ggplot2, ggExtra, dplyr, broom, shiny, shinyBS, shinycssloaders, shinydashboard, shinyjs, BiocStyle, rmarkdown, knitr, kableExtra, magick, VennDiagram, UpSetR, heatmaply, viridis, plotly, DT, gridExtra, methods, stats, SummarizedExperiment, BiocGenerics, S4Vectors, zip, reshape2

**Suggests** testthat

**SystemRequirements** To generate html reports pandoc (<http://pandoc.org/installing.html>) is required.

**biocViews** Software, TimeCourse, QualityControl, Visualization, GUI, PrincipalComponent

**URL** <https://github.com/matthewcarlucci/DiscoRhythm>

**BugReports** <https://github.com/matthewcarlucci/DiscoRhythm/issues>

**Config/pak/sysreqs** cmake libfontconfig1-dev libfreetype6-dev libfribidi-dev make libharfbuzz-dev libmagick++-dev gsfonts libicu-dev libpng-dev libuv1-dev libxml2-dev libssl-dev zlib1g-dev

**Repository** <https://bioc-release.r-universe.dev>

**Date/Publication** 2026-04-28 12:49:58 UTC

**RemoteUrl** <https://github.com/bioc/DiscoRhythm>

**RemoteRef** RELEASE\_3\_23

**RemoteSha** 390041f08f1f9f59d56aaea8d49131397091351c

## Contents

discoApp . . . . .	2
discoBatch . . . . .	3
discoCheckInput . . . . .	5
discoDesignSummary . . . . .	6
discoDFtoSE . . . . .	6
discoGetSimu . . . . .	7
discoODAexclusionMatrix . . . . .	8
discoODAIid2name . . . . .	9
discoODAs . . . . .	9
discoParseMeta . . . . .	11
discoPCA . . . . .	12
discoPeriodDetection . . . . .	13
discoQC . . . . .	14

<b>Index</b>	<b>16</b>
--------------	-----------

---

discoApp	<i>Launch the DiscoRhythm Shiny Application</i>
----------	---

---

### Description

This launches the web interface to DiscoRhythm containing all analysis tools. The vignette contains details on usage.

### Usage

```
discoApp(ncores = 1, port = 3838, local = TRUE)
```

### Arguments

ncores	numeric, number of cores to use for parallelized tasks. Currently, only used in oscillation detection function discoODAs.
port	numeric, port to run the shiny application on. Sets shiny.port option.
local	logical, set to FALSE for public server mode to reduce file size limits.

### Value

Nothing is returned by this function.

## Examples

```
## Not run:  
discoApp()  
  
## End(Not run)
```

---

discoBatch

*Core DiscoRhythm Workflow*

---

## Description

Execute the DiscoRhythm workflow with one command to obtain the results of oscillation detection (discoODAs) and optionally generate an html report with data visualizations from an Rmarkdown template. See the DiscoRhythm vignette for more details on the analysis procedures.

## Usage

```
discoBatch(  
  indata,  
  report = NULL,  
  outdata = TRUE,  
  ncores = 1,  
  timeType = "circular",  
  main_per = 24,  
  cor_threshold = 3,  
  cor_method = "pearson",  
  cor_threshType = "sd",  
  pca_threshold = 3,  
  pca_scale = TRUE,  
  pca_pcToCut = paste0("PC", seq_len(4)),  
  aov_method = "None",  
  aov_pcut = 0.05,  
  aov_Fcut = 0,  
  avg_method = "Median",  
  osc_method = NULL,  
  osc_period = 24  
)
```

## Arguments

indata	SummarizedExperiment or data.frame, see the vignette for the specific formats expected for each of these input types. discoParseMeta.
report	character, if !is.null(report) an html report with
outdata	logical, whether to return the final discoODAs (note if run with is.null(report) discoBatch will return nothing).

ncores	numeric, number of cores to use for parallelized tasks. Currently, only used in oscillation detection function discoODAs.
timeType	character, nature of the sample times provided (one of "circular" or "linear").
main_per	numeric, the length of the main hypothesized period (e.g. 24hr for circadian experiments). Used in discoPeriodDetection.
cor_threshold	numeric, threshold used in inter-sample correlation analysis for outlier detection. Either in units of correlation coefficient or standard deviations from the mean (see cor_threshType).
cor_method	character, which correlation method to use for outlier removal (see cor for more details).
cor_threshType	character, one of "sd" or "value" indicating whether cor_threshold should be set by absolute correlation coefficient or by standard deviations from the mean of all samples.
pca_threshold	numeric, the number of standard deviations to set as the threshold for outlier detection in PCA outlier removal.
pca_scale	logical, whether to scale the data prior to PCA.
pca_pcToCut	character, names of which PCs to use for outlier detection (e.g. "PC1", "PC2" etc.).
aov_method	character, method to use for ANOVA. One of: "Equal Variance", "Welch", or "None".
aov_pcut	numeric, p-value cutoff used to select rows with statistically significant signal-to-noise.
aov_Fcut	numeric, F-statistic cutoff used to select rows with high signal-to-noise based on magnitude.
avg_method	character, method for averaging technical replicates. One of: "Median", "Mean", "Random", or "None".
osc_method	character, vector of oscillation detection algorithms to apply to the data. Methods that are determined to be inappropriate for the experimental design (using the discoODAexclusionMatrix) will be ignored. If is.null(osc_method) all suitable methods will be executed.
osc_period	numeric, a fixed period to use for oscillation detection using all methods.

**Value**

returns the results of discoODAs

**See Also**

discoODAs, discoRepAnalysis, discoPeriodDetection, discoPCAoutliers, discoInterCorOutliers

**Examples**

```
indata <- discoGetSimu()

# Batch execute (on demo data) to generate a DiscoRhythm_report.html report.
```

```
# Returns the results of discoODAs
discoODAses <- discoBatch(indata,
  report="DiscoRhythm_report.html",
  osc_method="CS")
```

---

**discoCheckInput***Import Data for DiscoRhythm Analysis*

---

### Description

Performs various checks and cleaning operations on the input data.

### Usage

```
discoCheckInput(se, n_min_sample = 3)
```

### Arguments

<code>se</code>	SummarizedExperiment, the main data object used by DiscoRhythm expected to contain <code>se\$ID</code> , <code>se\$ReplicateID</code> , <code>se\$Time</code> sample metadata and non-null row-names. See the vignette for more details.
<code>n_min_sample</code>	numeric value specifying minimal number of samples needed to perform analysis.

### Details

Rows containing NA's or all constant values are removed. If matrix values are character it will be attempted to convert them to numeric. If input is not a matrix it will be converted using `as.matrix()`. User will be warned if row IDs contain duplicate entries.

### Value

SummarizedExperiment checked for errors and modified as needed

### Examples

```
se <- discoGetSimu(TRUE)
se_clean <- discoCheckInput(se)
```

---

discoDesignSummary      *Summarize the experimental design*

---

### Description

Using sample times and biological sample Ids, constructs a summary table of the number of total samples at each timepoint and additionally summarizes the number of replicates for each biological sample.

### Usage

```
discoDesignSummary(Metadata)
```

### Arguments

Metadata      data.frame of sample data, usually generated by using discoParseMeta on the column names of the Maindata data.frame. If is.null(Metadata) and Maindata is provided as input, Metadata will be generated from Maindata.

### Value

A table where the first row summarizes the number of datapoints for each timepoint and other cells indicate the number of technical replicates for a given biological sample.

### See Also

discoParseMeta

### Examples

```
# import example data
Metadata <- SummarizedExperiment::colData(discoGetSimu(TRUE))
# Summarize the experiment design
discoDesignSummary(Metadata)
```

---

discoDFtoSE      *Data formatting for DiscoRhythm*

---

### Description

Functions to import a data.frame (from the format expected by the web application discoApp()) as a SummarizedExperiment object or to export a SummarizedExperiment for use with the web application.

**Usage**

```
discoDFtoSE(Maindata, Metadata = NULL, shinySession = NULL)
```

```
discoSEtoDF(se)
```

**Arguments**

Maindata	data.frame with the first column containing row IDs and all subsequent columns containing experimental values. Columns should follow the expected naming format described in the vignette.
Metadata	data.frame of sample data, usually generated by using discoParseMeta on the column names of the Maindata data.frame. If is.null(Metadata) and Maindata is provided as input, Metadata will be generated from Maindata.
shinySession	shiny session object for use only by the DiscoRhythm shiny app discoApp() to update the axis labels using the time value prefix.
se	SummarizedExperiment, the main data object used by DiscoRhythm expected to contain se\$ID, se\$ReplicateID, se\$Time sample metadata and non-null row-names. See the vignette for more details.

**Value**

discoDFtoSE returns a SummarizedExperiment object with colData containing sample metadata.

discoSEtoDF returns a DiscoRhythm format data.frame.

**Examples**

```
df <- discoGetSimu()
se <- discoDFtoSE(df)

df <- discoSEtoDF(se)
```

---

discoGetSimu

*Read in the DiscoRhythm Simulated dataset*

---

**Description**

A convenience function to get the simulated circadian transcriptomic system data file used in DiscoRhythm for various demonstrations and tests.

**Usage**

```
discoGetSimu(as_se = FALSE)
```

**Arguments**

as_se	logical, indicates if example data should be returned as a SummarizedExperiment or data.frame.
-------	--

**Value**

The simulated demo dataset used in the DiscoRhythm web application as a data.frame or SummarizedExperiment.

**Examples**

```
indata <- discoGetSimu()
```

---

discoODAexclusionMatrix

*Algorithm Exclusion Matrix*

---

**Description**

A small matrix indicating which algorithms should be excluded given certain experimental designs and data types.

**Usage**

```
discoODAexclusionMatrix
```

**Format**

An object of class matrix (inherits from array) with 4 rows and 7 columns.

**Examples**

```
# Code used to generate discoODAexclusionMatrix

itemNames <- c(
  "missing_value",
  "with_bio_replicate",
  "non_integer_interval",
  "uneven_interval",
  "circular_t",
  "invalidPeriod",
  "invalidJTKperiod"
)

# Creating requirements matrix, first assuming all methods are valid
# Then applying exclusion criteria of MetaCycle plus CS criteria
mat <- matrix(TRUE, nrow = 4, ncol = length(itemNames))
rownames(mat) <- c("CS", "JTK", "LS", "ARS")
colnames(mat) <- itemNames

# Exclusion criteria from MetaCycle v1.1, i.e. can algorithm handle XXX
mat[c("ARS", "JTK"), c("non_integer_interval", "uneven_interval")] <- FALSE
```

```

mat["ARS", "with_bio_replicate"] <- FALSE
mat["ARS", "missing_value"] <- FALSE
mat["JTK", "invalidJTKperiod"] <- FALSE

# Additional exclusion criteria
mat["ARS", "circular_t"] <- FALSE
mat[c("CS", "JTK", "ARS", "LS"), "invalidPeriod"] <- FALSE

discoODAexclusionMatrix <- mat

```

---

discoODAid2name      *Mapping Identifiers to Full Names*

---

### Description

A small named vector mapping oscillation detection algorithm names to a convenient identifier.

### Usage

```
discoODAid2name
```

### Format

A named vector, length 4

**names(discoODAid2name)** Identifier

**as.vector(discoODAid2name)** Full names

---

discoODAs      *Execute Oscillation Detection Using DiscoRhythm*

---

### Description

Runs specified oscillation detection algorithms (ODAs) sequentially to obtain oscillation characteristics for each row of the input data.

### Usage

```
discoGetODAs(se, method = NULL, period, circular_t = FALSE)
```

```

discoODAs(
  se,
  period = 24,
  method = c("CS", "JTK", "LS", "ARS"),
  circular_t = FALSE,
  ncores = 1
)

```

**Arguments**

se	SummarizedExperiment, the main data object used by DiscoRhythm expected to contain se\$ID, se\$ReplicateID, se\$Time sample metadata and non-null row-names. See the vignette for more details.
method	character, short names of ODAs to use. If length>1 all input method names will be evaluated.
period	numeric, the hypothesized period to test for.
circular_t	logical, is time circular on some base-cycle (ex. time of day). See the DiscoRhythm vignette for details.
ncores	numeric, number of cores to parallelize with (applicable to JTK, ARSER and LS only). If 1, will execute in serial.

**Details**

There are currently 4 available algorithms for rhythm detection:

- CS = Cosinor (Cornelissen,G. 2014): a.k.a “Harmonic Regression” fits a sinusoid with a free phase parameter.
- LS = Lomb-Scargle (Glynn, 2006): an approach using spectral power density.
- ARS = ARSER (Yang, 2010): removes linear trends and performs the Cosinor test.
- JTK = JTK Cycle (Hughes, 2010): non-parametric test of rhythmicity robust to outliers.

LS, ARS, and JTK results come directly from MetaCycle meta2d() output using the specified fixed period. ARSmle is set to "nomle" and no method integration is used (see meta2d documentation for details).

CS is implemented directly in DiscoRhythm’s lmCSmat() as the single-component cosinor described in Cornelissen,G. (2014).

All q-values are calculated by performing p.adjust() on the resulting p-values with method="fdr".

Technical replicates are expected to be merged (likely by discoRepAnalysis) prior to usage of discoODAs.

The discoGetODAs function is called by discoODAs to determine if the selected methods may be used. If any methods are not valid, a warning will be thrown and only valid methods will be computed. discoGetODAs is not typically used directly, however, it may be called by the user to determine if the provided SummarizedExperiment is suitable for use with the specified methods.

**Value**

A named list of results where each element is a data.frame for the corresponding method with rownames corresponding to the feature identifiers and columns containing estimates for:

- acrophase
- amplitude
- p-value
- q-value

Additional columns relevant to each method will be present.

## References

- Yang R. and Su Z. (2010). Analyzing circadian expression data by harmonic regression based on autoregressive spectral estimation. *Bioinformatics*, **26(12)**, i168–i174.
- Hughes M. E., Hogenesch J. B. and Kornacker K. (2010). JTK\_CYCLE: an efficient nonparametric algorithm for detecting rhythmic components in genome-scale data sets. *Journal of Biological Rhythms*, **25(5)**, 372–380.
- Glynn E. F., Chen J. and Mushegian A. R. (2006). Detecting periodic patterns in unevenly spaced gene expression time series using Lomb-Scargle periodograms. *Bioinformatics*, **22(3)**, 310–316.
- Cornelissen, G. (2014) Cosinor-based rhythmometry. *Theor. Biol. Med. Model.*, **11**, 16.

## See Also

[lmCSmat meta2d](#)

## Examples

```
# Return valid ODAs for example dataset
discoGetODAs(discoGetSimu(as_se=TRUE),period=24)

# Import the simulated example dataset
se <- discoCheckInput(discoGetSimu(TRUE))

# Use discoRepAnalysis to average technical replicates
se_merged <- discoRepAnalysis(se,aov_pcut=1)$se

# Execute the Cosinor and JTK methods with a 24hr period
discoODAs <- discoODAs(se_merged,method=c("CS","JTK"))

# Get the index of rhythmic features detected by both methods at qvalue<0.05
idx <- which(discoODAs$CS$qvalue<0.05 & discoODAs$JTK$qvalue<0.05)

# Get the identifiers for common rhythmic features
rownames(se_merged)[idx]
```

---

discoParseMeta

*Generate Experiment Metadata*

---

## Description

Parses the sample metadata from a vector of sample names (often column names of a Maindata format data.frame).

## Usage

```
discoParseMeta(samplenames, shinySession = NULL)
```

**Arguments**

- `samplenames` character, a list of sample names following the DiscoRhythm naming convention (<prefix><Time>\_<UniqueID>\_<ReplicateID>).
- `shinySession` shiny session object for use only by the DiscoRhythm shiny app `discoApp()` to update the axis labels using the time value prefix.

**Details**

The regular expression used to obtain metadata is `"^([[:alpha:]]*)(\\-?[0-9]+[\\.]?[0-9]*)\\_? ([[:alnum:]]\\.)*\\_?([[:alnum:]]\\.)*$"`

Where each `()` will be used to construct the final metadata `data.frame`

**Value**

a `data.frame` containing 3 columns of metadata. `ID` = unique sample identity. `Time` = sample collection time. `ReplicateID` = Identifier where `Time + ReplicateID` indicates a biological sample ID.

**Examples**

```
discoParseMeta(c("CT24_AD_1", "CT24_AS_1", "CT24_AE_2", "CT24_AW_2",
"CT26_AB_1", "CT26_AC_1", "CT26_BB_2", "CT26_BC_2"))
```

---

discoPCA

*Perform PCA*

---

**Description**

Calculates PCA results from `prcomp` with error handling and outputs suitable for the DiscoRhythm workflow.

**Usage**

```
discoPCA(se, scale = TRUE, npcs = 10)
```

**Arguments**

- `se` SummarizedExperiment, the main data object used by DiscoRhythm expected to contain `se$ID`, `se$ReplicateID`, `se$Time` sample metadata and non-null row-names. See the vignette for more details.
- `scale` logical, whether or not to scale the data prior to PCA, see [prcomp](#) for more details.
- `npcs` numeric, maximum number of principal components to return.

**Value**

output from `prcomp` with an added table summary

**Examples**

```
se <- discoGetSimu(TRUE)
pca <- discoPCA(se)
```

---

`discoPeriodDetection` *Detect dataset-wide fits to multiple periodicities*

---

**Description**

Detect dataset-wide fits to multiple periodicities

**Usage**

```
discoPeriodDetection(
  se,
  timeType = c("linear", "circular"),
  main_per = 24,
  test_periods = NULL
)
```

**Arguments**

<code>se</code>	SummarizedExperiment, the main data object used by DiscoRhythm expected to contain <code>se\$ID</code> , <code>se\$ReplicateID</code> , <code>se\$Time</code> sample metadata and non-null row-names. See the vignette for more details.
<code>timeType</code>	character, time is either reported as "linear" or "circular" on some base-cycle (ex. time of day). This determines the periods that will be tested for.
<code>main_per</code>	numeric, if <code>timeType=="circular"</code> <code>main_per</code> indicates the period of the base-cycle where sampling times are derived.
<code>test_periods</code>	numeric, a vector of the periods to test. if <code>timeType=="linear"</code> and <code>length(test_periods)==2</code> it will be assumed to be a range of periods to test over.

**Value**

A data.frame of Rsquared values for each period, for each row of Maindata.

**Examples**

```
se <- discoGetSimu(TRUE)

# Detect periods
rsqs <- discoPeriodDetection(se)
```

**Description**

Functions for executing outlier detection and row filtering procedures prior to rhythmicity analysis.

**Usage**

```
discoPCAoutliers(se, threshold = 3, scale = TRUE, pcToCut = seq_len(4))
```

```
discoInterCorOutliers(
  se,
  cor_method = c("pearson", "kendall", "spearman"),
  threshold = 3,
  thresh_type = c("sd", "value")
)
```

```
discoRepAnalysis(
  se,
  aov_method = c("Equal Variance", "Welch", "None"),
  aov_pcut = 0.05,
  aov_Fcut = 0,
  avg_method = c("Median", "Mean", "Random", "None")
)
```

**Arguments**

se	SummarizedExperiment, the main data object used by DiscoRhythm expected to contain se\$ID, se\$ReplicateID, se\$Time sample metadata and non-null row-names. See the vignette for more details.
threshold	numeric, a threshold determining which samples are outliers (for discoInterCorOutliers, in units of thresh_type, for discoPCAoutliers in units of standard deviations).
scale	logical, whether or not to scale the data prior to PCA, see <a href="#">prcomp</a> for more details.
pcToCut	numeric, which PCs to use for outlier detection. It is recommended to select the first X PCs based on which PCs explain a significant amount of variance in the data.
cor_method	character, method of pairwise correlation (see <a href="#">cor</a> 's "method" argument for all options).
thresh_type	character indicating threshold type (either standard deviations below the mean, or an absolute correlation value). One of: "sd" or "value".
aov_method	character, method to use for ANOVA. One of: "Equal Variance", "Welch", or "None".

aov_pcut	numeric, p-value cutoff used to select rows with statistically significant signal-to-noise.
aov_Fcut	numeric, F-statistic cutoff used to select rows with high signal-to-noise based on magnitude.
avg_method	character, method for averaging technical replicates. One of: "Median", "Mean", "Random", or "None".

**Value**

list containing PCA results and the detected outliers

A list of 3 objects: 1) outliers - named logical indicating if the sample is an outlier 2) meanCor - mean of all pairwise correlations for a given sample 3) corMat - Matrix of all pairwise correlation values

**Examples**

```
se <- discoGetSimu(TRUE)
PCARes <- discoPCAoutliers(se)

CorRes <- discoInterCorOutliers(se)

ANOVARes <- discoRepAnalysis(se)
```

# Index

## \* datasets

- discoODAexclusionMatrix, 8
- discoODAid2name, 9

cor, [4](#), [14](#)

discoApp, [2](#)

discoBatch, [3](#)

discoCheckInput, [5](#)

discoDesignSummary, [6](#)

discoDFtoSE, [6](#)

discoGetODAs (discoODAs), [9](#)

discoGetSimu, [7](#)

discoInterCorOutliers (discoQC), [14](#)

discoODAexclusionMatrix, [8](#)

discoODAid2name, [9](#)

discoODAs, [9](#)

discoParseMeta, [11](#)

discoPCA, [12](#)

discoPCAoutliers (discoQC), [14](#)

discoPeriodDetection, [13](#)

discoQC, [14](#)

discoRepAnalysis (discoQC), [14](#)

discoSEtoDF (discoDFtoSE), [6](#)

lmCSmat, [11](#)

meta2d, [11](#)

prcomp, [12–14](#)